#### DENTSLEEVE MANOMETRIC CATHETERS

#### INSTRUCTIONS FOR USE

#### Measure to be Taken Prior to Intubation

- Check that the design of the manometric catheter is specified for the intended measurements. Catheters are designed for recording from specific regions.
- Only use catheters that have been cleaned and high-level disinfected (or sterilized) according to instructions below.
- Inspect the catheter for any evidence of mechanical damage or imperfections. **Do not use if damaged.**
- Water fill channels that will not be used for measurements and seal catheter luer connectors with the luer plugs provided.

#### Other equipment and Expertise Necessary for Measurements

- This catheter must be perfused with a manometric perfusion pump, at the recommended flow rate for achievement of technically valid measurements.
- Pressure amplification and recording equipment are also essential.
- This catheter must only be used under medical supervision by personnel who have been adequately trained in safe gastro-intestinal intubation, and in the performance of perfused manometric measurements.

#### Passage and Positioning of Catheter

- Dentsleeve upper gastro-intestinal manometric catheters are designed primarily for trans-nasal passage, after application of local anaesthetic.
- In the case of pharyngeal and oesophageal manometric catheters, correct positioning of the catheter is determined by evaluation of the pressure tracings.
- Special techniques are needed for reliable passage of manometric catheters across the pylorus and along the small intestine; verification of catheter position also requires special measures.

#### Positioning of Sleeve Catheters (if applicable)

- The Sleeve should be positioned so that it remains within the sphincter during its normal movement.
- After deeper insertion of the Sleeve, pull the catheter back into the sphincter, observing tracings from the Sleeve and side holes at either end of the Sleeve.
- Recording of a sphincteric pattern of pressure from a side hole at either end of the Sleeve indicates incorrect placement of the catheter relative to the sphincter.

#### Cleaning and Sterilisation

This catheter is a re-useable device which has been validated as being able to withstand 50 cycles of steam autoclaving.

- Immediately after extubation, immerse the catheter in a bowl of warm, mild detergent solution and wipe externally

145 Traders Blvd. East, Unit 33-34, Mississauga, Ontario, Canada L4Z 3L3
Tel: (905) 712-8828 Fax: (905) 890-3523 Toll Free: (866) 751-6345
www.dentsleeve.com Email: mail@dentsleeve.com

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several times. Do not allow the uncleaned catheter to dry out.

- Luer connector plugs must be removed for cleaning and autoclaved with the catheter.
- Flush the detergent solution gently through each channel with a 20ml syringe – Do not over-pressurise.
- Rinse the outside of the catheter in a bowl of clean water.
- Place the catheter on a clean towel, and cover it with a fold of the towel.
- Flush each channel first with water, then air – do not over-pressurise.
- Package the catheter for autoclaving, and autoclave without delay at 134°C for five minutes at 206 kPa or 30 psi, with a total cycle time of 30 minutes to allow for warm up and cool down time.

**Intended Use**

- Dentsleeve Manometric Catheters have been designed and manufactured for measurement of intraluminal gastro-intestinal pressures in both adults and children when used with a manometric perfusion pump.
- This manometric catheter is to be used only under the supervision of a physician who has received adequate training in gastro-intestinal manometry.

**Indications for Use**

- Use of Dentsleeve Manometric Catheters is indicated when measurements of gastro-intestinal tract pressures are judged to be useful for determining management of patients with proven or suspected gastro-intestinal motor disorders.

**Contraindications for Use**

- Where there is a significantly increased risk of pulmonary aspiration associated with intubation which can not be adequately reduced by practical measures such as a longer than usual period of fasting prior to intubation.
- The presence of any stenosis, structural deformity or disease that results in an unacceptably increased rick of aspiration, perforation or any other adverse consequence arising from use of this assembly.

**WARNING**

- This devise is supplied non-sterile.
- Incorrect manometric perfusion pump operation can cause delivery of large volumes of gas through a manometric catheter during measurements.
- The use of "off-the-shelf" components from other manufacturers is not recommended. Only Dentsleeve components have been tested and validated to work safely and accurately with this catheter.
- Safe use of this device requires supervision by a physician trained in gastro-intestinal manometry.

**USA only**

***Caution: US Federal law restricts this device to sale by or on the order of a physician.***